

— Soft laser expertise

# premio 30

Laser duo



Acupuncture  
Ear Acupuncture

Sedatelec


# premio 30

Laser duo

Painless stimulation of body and auricular points



**4 treatment  
modes**



**2 powers,  
adaptable for  
superficial or  
deeper points**



**10 hours  
working time**  
Rechargeable  
battery

# Between tradition and modernity: **safe and effective**

## Acupuncture

### Tonification

Controlled, rapid  
tonification  
of acupuncture points.



### Dispersion

In a few seconds,  
**the Qi will be disperse**  
and you will remove  
energy stasis.



### Harmonisation

This treatment **rebalances  
the point and the circulation  
in the meridian**  
by alternating impulse rates.



## Ear Acupuncture

### Frequency stimulation

This selection  
will guarantee stimulation  
of **auricular** points.




# Innovation and effectiveness of frequencies

## New impulse emission dynamics

### Wobulation (100-10 000Hz)


#### Exponential acceleration

Tonification occurs as a result of an **acceleration** in **impulse emission rate**, which has an energising effect.




#### Exponential deceleration

Dispersion is obtained by **slowing** the **impulse emission rate**.



#### Alternating acceleration and deceleration

Alternating **acceleration** and **deceleration** of the emission rate leads to harmonisation of the point.



## Scanning

### According to the 7 Nogier\* frequencies

The message is contained in the frequency for this ear acupuncture programme. It provides a successive emission by **scanning all the Nogier\* frequencies** (from A to G, 73 - 4 672 Hz).


*\* The application of Nogier frequencies with this device has not been authorised by Health Canada for lack of sufficient clinical studies.*




# Renowned experience in Pulsed Laser

## Wavelength

The **905nm infra-red emission** passes through the skin barrier and reaches the deep cells, entering entirely safely without being absorbed by skin water, haemoglobin or melanin.



## Peak value



The laser beam is emitted as short impulses of high peak power. This method of emission gives the laser beam the **necessary power to enter tissues** entirely safely with no destructive heat effect, carrying the useful information to the site point.

## Adjustable power

15 or 40 Watts power is adjustable depending on the desired depths

- **Superficial points:** surface points (face, ears)
- **Deep points:** other body points

## Application time

**Short treatment of 30 seconds** for each point. This application time provides the optimal energy per point taking account of the physiological body reaction time.

**Sedatelec +:**  
**frequencies wobulation for acupuncture**  
The Sedatelec signature:  
**Nogier\* frequencies for ear acupuncture**



# premio 30

Laser duo

**Effective, light  
and ergonomic**

**A wireless laser,  
ideal for the  
consulting room  
and on the move**



*Invisible laser emission. Do not point the laser beam towards eyes. 3B classification according to IC60825-1*

*Caution : Federal law restricts this device to sale by or on the order of a physician.*

Your distributor

## **Eastern Currents**

9109 Shaughnessy St

Vancouver, BC V6P 6R9

**[www.easterncurrents.ca](http://www.easterncurrents.ca)**

**1-800-667-6866**



Made in France

**Sedatelec** • Chemin des Mûriers • F-69540 Irigny • Lyon • France • Tél. +33 (0) 472 663 322  
Fax + 33 (0) 478 508 903 • [sedatelec@sedatelec.com](mailto:sedatelec@sedatelec.com) • **[www.sedatelec.com](http://www.sedatelec.com)**

